



Midmark 223 Exam Table- Fully Refurbished

Getting on an exam table can be an anxious, insecure moment, especially for elderly and special needs patients. Ritter 223 Barrier-Free® Tables are designed to lower to 18 inches from the floor, making it easier for patients to get on the table and eliminating excess strain on their backs as well as the backs of physicians and staff. This low height means patients can get on the table with no or minimal staff assistance, helping them maintain their dignity.

Features:

- Power height: 18" to 37" seat adjustment
- Power back with spring air assist
- 400 lb patient weight capacity
- Foot control
- Stirrups
- Pelvic tilt
- Drawer heater
- Paper roll holder (hidden, holds 21" x 3.5" paper roll)
- Electric receptacle (115 VAC 5 amps max) duplex, hospital grade located on the left side of table
- Retractable foot holders
- Ritter vinyl seamless unholstery
- Foot control only

Specs:

- Length with footrest extended: 72" (182.9 cm)
- Height: Minimum 18" (45.7 cm) Maximum 37" (94.0 cm)
- Patient Load Rating 400 lb (181.4 kg)
- Paper roll holder (hidden): Holds 21" x 3.5" paper roll (53.3 x 8.9 cm)
- Electrical receptacle (115 VAC 5 Amps Max.) Duplex, hospital grade, electrical receptacle on left side



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.