



Midmark 244 Bariatric Power Treatment Table - Fully Refurbished

Midmark 244 Bariatric Treatment Table is successfully addressing special needs of certain patient groups, ensuring their safety, comfort and confidence. It is designed to aim most demanding medical environments minimizing the pressure on staff and delivering a demanded functionality to a practitioner. Experience the benefits of an improved accessibility, vast multipositioning capabilities and ease of operation. The 850 pound weight capacity, 18 inch low height, 34 inch high height and 32 inch wide seat make patient transfer a very simple routine, requiring minimal assistance from you and your staff, as well as making this the perfect table for wound care procedures.

Features:

- Power Height, Back and Tilt
- Powered multi-positioning capabilities
- Reverse Trendelenburg
- Ergonomic, intuitive standard hand control in an easy-to-find location
- Back support angle: 0° to 65°
- Seat tilt angle: -15° to 15°
- Paper roll holder capacity: 18", 21", or 24"
- Ritter 244 Bariatric Power Treatment Tables meet UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards

Specs:

- 850 Pound Patient Weight Capacity
- 18 inch to 34 inch Seat Height Adjustment
- 32 inch Wide Seat
- Head and foot area: 30" W (76.2 cm)
- Overall length: 80" (203.2 cm)



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.