



Midmark 510 Stretcher- Fully Refurbished

Midmark Stretchers are one of the most reliable and well known brands of stretchers available. Midmark offers stretchers for many areas of your medical facility, including transport, PACU, trauma eye/ophthalmic, pediatric, and many more. Midmark stretchers have a single pedestal hydraulic design. Most Midmark stretchers have a central brake and steer function, gas-lift assisted fowler (backrest), trendelenburg / reverse trendelenburg and hydraulic height adjustment. Fold under side rails provide maximum access to the patient during treatment. While new stretchers are not in production by Midmark, used Midmark stretchers and refurbished Midmark stretchers can be a great option to purchase a quality Midmark stretcher at a very affordable price.

Features:

- Basket for Patient's belongings
- IV Pole
- Pneumatic Fowler Position
- Trendelenburg: 12°
- Reverse Trendelenburg: 7°
- Fowler: 90°
- Back raises and lowers (manual – hydraulic assist), whole table tilts.
- Easy to raise and lower hospital side rails.
- Easy to use "One lock – all lock" wheels.
- 30" W Total
- 25" Wide Patient Surface
- Unit Raises and lowers in trendelenberg only (not reverse trend)
- Back Raises and Lowers – Hydraulic
- Weight: 107 kg
- Central brakes
- Conductive Wheels
- 6 sockets
- Locking pins
- Steer and Side Rails



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.