



Midmark 625 Power Exam Table with Scale- Fully Refurbished

The Midmark 625 with IQscale™ features a fully integrated and accurate digital scale, which weighs patients quickly, efficiently and discreetly. All-in-one handheld control does it all. This simple handheld device controls the entire table, including seat height and seatback angle adjustment, scale functions and BMI calculations. Since patients remain positioned on the exam table, you'll streamline your workflow, save valuable time, and be able to spend more time caring for your patients. The Midmark 625 with IQscale™ allows physicians to care for a growing number of patients without worrying about equipment limitations.

Features:

- Power height and back
- 650-pound patient weight capacity
- 18-inch to 37-inch seat height adjustment
- Unique Quick Exam® and Motion Profiling™ control features for more efficient patient positioning
- Optional IQscale™ and IQHub™ offers a digital scale built into the table and easy access receptacles and data ports
- Powered Multi-Positioning
- Standard hand and foot control make patient positioning quick, effective and always within reach.



Specs:

- Height: Minimum 18" (46 cm) Maximum 37" (94 cm)
- Premium Upholstery Dimensions (top): 28" W x 59.5" L (71 cm x 151 cm)
- Length with footrest extended: 75" (190 cm)
- Stirrups are standard
- Pelvic Tilt: 7°
- Paper Roll holder: holds 21" x 3.5" diameter paper roll (53 cm x 9 cm)
- Electrical Receptacle: 115 VAC, 5 A Max.

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.