



Midmark 647 Power Podiatry Chair - Fully Refurbished

The 647's new, ergonomic foot section design gives physicians what they need at the point of care to perform even the lengthiest procedures. Touch controls and sealed foot control are standard on this model, as are covered glides on the foot section and debris tray. The easy-to-reach, release handles allow for one-handed, smooth, quiet operation. Integrated touch controls means there is no hand control to get knocked to the ground, eliminating costly breakage and allowing unrestricted patient access. You'll be surprised at how easy it is to keep your chair clean and operating smoothly.

Features:

- Power Tilt and back
- 450 lb. weight capacity
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Armrests standard
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards

Specs:

- Base height: Minimum 19"; Maximum 40" NOTE: Rotation adds 3.5" to minimum and maximum seat heights
- Foot section: 26"W x 6"H
- Foot extension: 9"
- Debris tray: 15.1"W x 10.6"H x 1.5"D
- Back support angle range: 0° to 80°
- Seat tilt angle range: 0° to 30°
- Seat width with arms installed – 37.5"
- Includes upholstery color of choice
- Table weight: 323 lbs



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.