



## Midmark Ritter M9, Old Style - Fully Refurbished

The M9 autoclave couldn't be easier and is time tested for reliability. You simply load it – set it – and forget it. Select your cycle and the time is set for you. Press "Start" and the M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam, giving you fast drying time for your instruments.

### Features:

- Fills with correct amount of water when you hit the start button.
- Door opens automatically to dissipate steam and give you the fastest drying time.
- Built-in safety features such as:
  - Overheat protection
  - Pressure relief valve
  - Flashing "Error" light.
- Pre-Programmed Controls for convenience
- Exclusive automatic-opening door for ease of use and a quick drying time
- Steam-Flush Pressure Pulse Air Removal System ensures proper steam penetration



### Specs:

- Overall length w/plug: 20.4" (51.8 cm)
- Overall width: 15.3" (38.9 cm)
- Overall height w/printer: 15.8" (40.1 cm)
- Minimum Countertop area: 15.3" x 17.9" (38.9 cm x 45.4 cm)
- Chamber:
  - 9" diameter (22.9 cm)
  - 15" depth (38.1 cm)
  - 3.5 gal. usable volume (13.4 liters)

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.