



auxoMEDICAL

NUVO Vu LED Surgical Light - New

When you're performing surgery, it's critical that you can see what you're doing down to the most minute detail. That's why the Nuvo Vu LED Surgical Light should be at the top of the list for your hospital operating rooms. The Nuvo Vu LED Surgical Light is brighter, whiter, cooler and more capable in the operating room than any previous technology. So you get the most accurate colors with unparalleled shadow control — and virtually no heat emission. All this for a value that's hard to match.

Features:

- Adjustable color temperature (3,800° to 4,800° K)
- Brilliant 160,000 lux output
- 50,000 hour-rated LED life
- Cross Focus Technology provides superior shadow control
- Near zero radiated heat
- Complete range of mounting options
- "Green Mode" for low-light, minimally- invasive surgery

Options:

- Solo Ceiling Mount
- Duo Ceiling Mount
- Trio Ceiling Mount

Optional "In-Light" HD camera

- 1080i high-definition
- HD-SDI video output
- 20x optical zoom/12x digital zoom
- Auto and manual focus
- Sterilizable handle - standard
- Operates from remote Wall Control



Why purchase equipment from Auxo Medical?

At Auxo Medical, we don't take any decision lightly. Whether we focus on quality of new or refurbished medical equipment, repair savings, or improved operational cost-efficiency—our eye is always on the smallest details that can affect the level of care provided and your budget. Our consultants pride themselves in helping you plan and develop the most appropriate facility for your changing needs.