



Passport 2 Anesthesia Monitor/5 Gas Module - Fully Refurbished

The Passport 2 makes monitoring simple. Rather than rummaging through layers of menus, there's a Navigator™ knob. For frequently used functions, there are quick-action keys. Three waveforms will automatically reconfigure to six when needed. Colorful numerics automatically re-adjust to maximum size for optimal display visibility. Be amazed at what a difference a truly intelligent interface will make. Who would have thought it possible? A monitor that actually provides more while requiring less.

Features:

- Multi-Parameter
- Measures ECG, NIBP, SpO2, Temp
- Gas Module for Co2 and AA
- Auto Agent Identification
- 3 or 5 Lead ECG
- Nellcor or Masimo SpO2
- Onboard Thermal Printer
- Fully User Configurable Display
- 4 Auto-Set Alarm Functions
- EMR Compatible
- Battery Operated w/ Built-In Charger
- Large, Bright 10.4" Color Display
- Lightweight Design
- Carry Handle
- 5 agent analysis

Specs:

- Dimensions: ~14"H x 12"W x 11"D
- Weight: ~24 Lbs
- Battery Voltage: 12V (2)
- Power Requirement: 120v
- Battery Run Time: Approx 2 Hours,
- Recorder Paper: 80mm, Thermal
- Battery: Rechargeable, Sealed Lead



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.