



# auxoMEDICAL

## Passport 2 with CO2 card - Fully Refurbished

The Datascope Passport 2 makes monitoring simple. Rather than rummaging through layers of menus, there's a Navigator knob. For frequently used functions, there are quick-action keys. Three waveforms will automatically reconfigure to six when needed. Colorful numerics automatically re-adjust to maximum size for optimal display visibility.

### Features:

- Dedicated function keys provide one-touch access to all vital functions
- Colorful waveforms with automatic display configuration simplify use and enhance display visibility
- Six waveforms
- Auto-adjustable large numerics
- 3 or 5-lead ECG
- Lead-selectable respiration
- Non-invasive blood pressure
- Temperature
- Masimo SET motion tolerant SpO2
- Optional features include: arrhythmia analysis, 3-lead ST analysis with trends and alarms, Microstream CO2, Nellcor OxiMax SpO2, two invasive blood pressures, and a dual trace recorder
- View 12 ECG analysis module offers diagnostic quality 12-lead ECG interpretation with ST and arrhythmia analysis in a compact, easy-to-use module
- Built-in PCMCIA ports allow you to upgrade software, enhance memory, and transfer data with minimal effort
- Integrated power supply with lithium-ion battery technology provides uninterrupted patient care during transport.
- CO2



\*Mounting Equipment Sold Separately

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.