



auxoMEDICAL

Philip SureSign VM6 - Fully Refurbished

The Philips SureSigns VM6 is designed to make patient care efficient and friendly. The VM6 is designed for easy mobility, with its easy to carry handle light weight, 7.1 lbs, and 4 hour battery life. Built with Philips technology the VM6 puts all patient vitals and monitoring at your finger tips. The VM6 performs pulse oximetry, non-invasive blood pressure, respiration, basic arrhythmia, and functions as a 3 or 5 lead ECG. There is also a USB port for data export and software upgrades. It has a lithium-ion battery as an internal power supply, and an AC power cord for operation and charging. Flexible and compact, the VM6 is ideal for continuous monitoring multiple vitals for adults, pediatrics, and neonates.

Features:

- Bright, color backlit 8.4" LED screen
- Audible & visual alarms with Nurse Call Output
- 4 display modes: Numeric, Numeric w/ Waveform, Tabular Trend, Graphical Trend
- Front panel icon shows power/battery status and battery charging
- Storage room for up to 96 hours of measurements
- Uses the same supplies as all other Philips monitors
- Intuitive Navigation; making frequently used test easily available
- Optional Recorder
- ECG, NIBP, SpO2, TEMP & printer
- Roll stand or wall mount option
- Dimensions: 26 cm W x 22 cm H x 14.5 cm D (10.2" W x 8.6" H x 5.7" D)
- Screen type: 21.3 cm (8.4") SVGA TFT-AM LCD display



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.