



Philips HeartStart FR-2 AED - Fully Refurbished

Heartstart FR2 provides emergency first responders with an effective means to treat victims in Sudden Cardiac Arrest (SCA). This lightweight, rugged and easy-to-use Automated External Defibrillator (AED) allows first responders to deliver lifesaving therapy as soon as they arrive on the scene. Its patented SMART Biphasic waveform is the most patient-friendly technology available today.

Features:

- Long life battery: Minimum 300 shocks when new. Shelf life before installation and standby life after installation typically 5 years each.
- Rechargeable training battery: A cost efficient way to utilize the device for training
- Backlit display
- Improves the quality of ECG and text prompts in conditions of direct sunlight and darkness
- Selftest and status indicator Monitors readiness and clearly displays whether or not the device is operational
- Easy and flexible configuration
- Allows extensive configuration possibilities, including treatment protocol in accordance with new ERC or AHA guidelines
- Manual override
- Adjust the contrast of the screen display and control special functions
- Data card tray
- Data card port
- Infrared (IR) communications port
- ECG display
- Case
- Pads



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.