



## Philips MRX Basic With SPO2 - Fully Refurbished

The Philips Heartstart MRx Defibrillator is the latest technology in biphasic defibrillation from Philips, a leader in patient monitoring. Philips uses the patented SMART biphasic waveform technology in this portable unit. The MRx is suitable for EMS and hospital use, as well as small surgery centers. Several features include the LONGEST BATTERY LIFE available in any comparable monitor, as well as the LARGEST COLOR DISPLAY and the fastest shock response time in the industry

### Features:

- Manual defibrillation
- AED mode
- ECG arrhythmia detection
- Automated testing
- Charge button charges the defibrillator in <5 seconds
- Large color display shows 4 waveforms and numerics, or view all 12 leads at once with the 12-lead acquisition option.
- Normal or high-contrast view for easy viewing in sunlight conditions.
- 10 hours of monitoring with two fully charged batteries.
- Ready-for-use indicator flashes to signal that ample power exists to monitor and deliver a shock.
- Automated self-tests that run hourly, daily and weekly
- Color-coded monitoring ports for easy identification
- Integration and upgrades made easy
- Storage capacity for up to 50 12-lead EKG recordings
- Temperature capability
- Ventilation monitoring during a code
- Analysis on screen
- CPR measurement and feedback
- SPO2



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.