



Physio-Control LIFEPAK 15 - Fully Refurbished

Physio-Control® defibrillators have set the standard for over 55 years, and the latest version of the LIFEPAK® 15 monitor raises the bar. The LIFEPAK® 15 device offers a sophisticated clinical technologies and operational effectiveness with a rich array of features – like the most powerful escalating energy available (up to 360J), advanced monitoring parameters, and a completely upgradable platform. At the same time, it's tough enough to stand up to your most challenging environments.

Features:

- Windows based organization for easy navigation
- A tiered view of all available setup options with defaults highlighted in bold
- Password access to prevent unauthorized changes
- Electronic and printed report options
- Additional language and device support files downloaded as needed
- ADAPTIV™ Biphasic Technology
- Manual & AED External Defibrillation
- 8.4" Color LCD Screen
- CPR Metronome
- Noninvasive Pacing (Demand & Non-Demand)
- 3-Lead ECG Monitoring
- 100mm Printer
- New monitoring parameters
- Advanced support for treating STEMI patients
- Most potent escalating energy available
- Proven CPR guidance
- Upgradable platform
- Latest Lithium-ion battery technology
- Toughened inside and out
- Dimensions: 12.5"(H) x 15.8"(W) x 9.1"(D)
- Weight: Basic monitor/defibrillator with new roll paper and two batteries installed: 17.5 lb



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.