



Puritan Bennett 7200AE Ventilator - Fully Refurbished

The Puritan Bennett 7200 respiratory ventilator has an advanced micro-processor design that offers CMV, assist/control, SIMV, CPAP, PEEP, and sigh controls. This critical care ventilator has a built-in nebulizer, air/oxygen blender, and variable pressure waveforms. Data is clear on the digital display and easy to input with the keyboard control.

Features:

- State-of-the-art gas delivery system.
- Demonstrated track record of reliability
- Full-featured microprocessor-controlled ventilator with 7202 display screen.
- Gas delivery features - Flow-by®, Pressure Support, Pressure Control Ventilation
- Patient monitoring features - Graphics 2.0, Respiratory Mechanics, Auto-PEEP
- Communications features – RS-232, Digital Communications Interface (DCI 2.0)
- Optional 7250⁺ Metabolic Monitor
- Integrated safety system
- Alarm Control Functions: Alarm silence, Alarm reset
- Parameters Display:
 - Airway pressure
 - Exhaled volume
 - Type of breath
 - Mean airway pressure and peak airway pressure
 - PEEP/CPAP pressure and plateau pressure
 - Respiratory rate
 - I:E ratio
 - Tidal volume and minute volume
- Dimensions: 16.5" (H) x 22" (W) x 22.5 (D)
- Weight: 50.8 kg (50.8 kg)
- Respiratory Rate: 0.5 to 70 bpm;
- Tidal Volume: .1 to 2.5 liters



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.