



Respironics PLV-102 Ventilator- Fully Refurbished

Respironics LifeCare PLV-102 Ventilator (with internal oxygen blender) – The LIFECARE PLV-102 is a microprocessor controlled, electrically powered portable volume ventilator that may also be pressure limited. Specifically designed for long-term use in the home, it is equally suitable for transport, short-term and institutional use. The ventilator is small and lightweight and will operate using a 12 volt DC external battery as well as 120 volt AC power. Reliability, versatility, and ease of operation are all hallmarks of this machine, making it the ideal choice for this setting.

Features:

- Can ventilate a wide variety of patients
- Low flow rates: 10L/min
- Low tidal volume: 50 ml
- Peak flow ability: 120 L/min
- Digital Displays:
 - Tidal Volume
 - I:E Ratio
 - Respiratory Rate
 - Inspiratory Flow Rate
- Alarm Interface Capabilities:
 - Remote Alarm
 - Pressure Alarm
 - Hospital/Nurse Call System
- Three Modes Of Operation:
 - Control
 - Assist-Control
 - Spontaneous Intermittent Ventilation (SIMV)
- Multiple Power Sources:
 - 120 VAC/240 VAC
 - 12 VDC external capabilities (up to 20 hours)
 - Internal battery (1 hour fully charged)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.