



Skytron 3100 - Fully Refurbished

The Skytron 3100 General Purpose Surgical Table has been optimized for all procedures requiring optimal C-arm access. The 23 ½" Sliding Top is a key feature of this table making it well suited to lower abdominal procedures such as GYN, Urology, Hip Revisions and Ureteroscopies. Despite the longitudinal travel of 23 ½" this table is still capable of supporting a full 500 lb patient regardless of the position required for surgery. The capabilities of the table have been further enhanced by the 7 inch lift provided by the Kidney Bridge.

Features:

- Radiographic top
- Precision locator system
- C-arm compatible
- Removable leg extension
- Remote control positioning
- Manual brake release
- Returns to level
- Standard side rails
- Cassette Tunnel

Specs:

Table Height: Raised 40" Lowered 29"

Trendelenberg and Reverse Trendelenberg: 20°

Lateral Tilt: ±20°

Kidney Bridge: 7"

Elevation Top Slide: 23 1/2"

Weight Limit: 500 lbs

Length: 77 "

Width: 19 5/8"

Width w/Rail: 21 7/8"

Weight of Table: 774 lbs



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.