



Skytron 3500 Elite - Fully Refurbished

The Skytron 3500 General Purpose Surgery Table takes patient positioning and imaging access to whole new level. The 21 inch sliding top combined with full patient articulation makes this the table of choice for virtually every major general surgical procedure. The low height has also been enhanced, making it better suited to procedures with the patient in the seated position as well as bariatric patients. This table has the ability to achieve extreme positions while maintaining patient safety. A stable footprint and rigid surface provides assurance that the patient position will be maintained as desired. The Skytron 3500 Elite features an incredible 21" top slide, 500 lb lift/articulation, return-to-level and return-to-center. An ultra modern table, able to satisfy all of your C-Arm and general surgery needs. Good for use in urology, neurosurgery, orthopedic, cysto, vascular, E.N.T. and open heart surgery.

Features:

- 500 pounds weight capacity
- Provides full body imaging
- Instant response table movements
- Easy-to-use pendant control

Specs:

- Trendelenburg: 25°/20°
- Reverse Trendelenburg: 25°/20°
- Lateral Tilt: 20°/20°
- Elevation: 26" – 45"
- Back Up/Down: 90°/40°
- Leg Up/Down: 0°/100°
- Leg: Removable
- Top Slide Head/Feet: 7½"/13½"
- Top Slide Total: 21"
- Return to Level: Returns all functions to level
- Kidney Elevator Up/Down: 5½"/0"
- Flex/Reflex: 120°/90°
- Emergency Brake Release: Fully Mechanical



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.