



Skytron 3501B - Fully Refurbished

Skytron 3501B EZ Slide is a universally versatile surgical table providing excellent imaging flexibility, clinical versatility and cost effectiveness. The EZ Slide 3501B Specialty table provides an impressive range of full-body imaging and articulating functions without the need to reposition the patient or resort to multiple support accessories. It features a 45-inch maximum and 26-inch minimum table height with 21-inch top side travel.

Features:

- 700 pound lift capacity.
- 600 pound articulate capacity.
- 21 inch longitudinal top slide.
- 7-1/2 inch slide to head.
- 13-1/2 inch slide to foot.
- 45 inch maximum table height.
- 26 inch minimum table height.
- Power kidney bridge.
- One touch (v2) current model backlit pendant control.



Specs:

- Length w /o head rest: 69 in. (1750mm)
- Length with head rest: 82 in. (2080mm)
- Width (back / seat plate): 19.75 in. (500mm)
- Width with side rails (back / seat plate): 21.75 in. (550mm)
- Weight: 827 lbs. (376kg)
- Permissible Load: 700 lbs. (320kg) Lift / 600 lbs. (275kg) Articulation
- Lateral tilt (left/right): 20°
- Trendelenburg: 25°
- Reverse Trendelenburg: 20°
- Power Kidney Bridge: 6 in. Lift 150mm Lift
- Back up/down: 90°/40°
- Leg up/down: 2°/98°

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.