



Skytron 6001 Elite with Battery Back-up - Fully Refurbished

The Skytron 6001 Elite Surgical Table is a remote-controlled operating table that allows for multi-positions of the patient and C-arm. Due to the power functions, this table is suitable for a wide range of procedures. Features include a 30 degree lateral tilt, trendelenberg and reverse trendelenberg, 180 degree top rotation, electro-hydraulic floor locks and a standard 500 lb weight capacity.

Features:

- Battery Backup
- Electro-hydraulic floor locks
- C-Arm compatible radiolucent top
- Remote control
- Manual Brake Release
- Kidney lift

Specs:

- Minimum Height: 28"
- Maximum Height: 45"
- Table Width: 20"
- Table Length: 75"
- Lateral Tilt: 30°
- Trendelenburg: 30°
- Top Rotation: 180°
- Weight Capacity: 500lbs



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.