



auxoMEDICAL

Skytron 6500 - Fully Refurbished

The Skytron® 6500 Surgical Table is a general surgery table filled with features. The 500 lb. (Elite) and 850 lb. (Hercules) weight capacity and the C-Arm compatibility allow this table to be put into any operating room. This table can also adjust to the trendelenburg, back, side, foot, and flex positions.

Features:

- 500 lb capacity
- C-Arm compatible
- X-ray top set
- 180 degree table top rotation
- Removable powered leg section
- Accommodates taller patients with the adjustable head section
- Use the reversible head section as a foot extension while transferring patients
- Sturdy, easy-to-clean stainless steel column and base
- Brakes are completely mechanical and provide dependable performance, even when there is no power to the machine
- Operates on line or battery power
- Push button break system
- Large perineal cut-out
- Kidney bridge with 5" of additional lift
- Height Range: 47" -26"
- Reflex: 90°
- Flex: 70°
- Foot / Leg Range: 105°- 20°
- Back range 40°-90°
- Trendelenburg and Reverse Trendelenburg: 30°
- Tilt Right/Left: 30°



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.