



Skytron 6700 - Fully Refurbished

The Skytron 6700 general surgical table represents the new standard in 21st Century Surgical Table design and performance, with unmatched imaging flexibility and even greater lift and articulation power. With its 1,200 lb patient weight capacity and the ability to be widened up to 16", it's a very popular operating table for bariatric surgery.

Features:

- Optional Infrared (IR) Remote Wireless Control
- 1,200 lb. lift & 1000 lb. articulate capacity
- 210¼ Top Rotation
- 30¼ patient articulation
- Return-to-level
- Auto Beach Chair Position
- Power kidney bridge
- Removable back & leg sections
- Optional Carbon Fiber back/leg sections
- Durable one-touch, backlit pendant control
- Battery power



Specs:

- 44" max./23" min. table height
- Lateral tilt of 30°
- Trendelenburg and reverse Trendelenburg: 30°
- Kidney bridge elevation: 5-3/4"
- Back up: 90°
- Back down: 40°
- Leg section down: 100°
- Leg section up: 20°
- top rotation: 210°
- 40° flex and 90° reflex function control button.

Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.