



Skytron 6701 - Fully Refurbished

From ease of set-up to the reduced need for multiple support accessories, the Hercules 6701 is the most powerful, versatile, and image-friendly surgical table available. The Hercules 6701 surgical table represents the new standard in 21st Century Surgical Table design and performance, with unmatched imaging flexibility and even greater lift and articulation power. The Hercules 6701 is a surgery table that can be used for general surgery or for bariatric patients.

Features:

- Capacity lift and articulate 1,200/1,000 pounds
- Height 23 low / 41 high inches, (6700B is 44 inches high).
- Auto Beach chair
- 35 degree rotation for ENT
- 210° table top rotation
- Lateral tilt right/left 30°
- Flex and reflex 40°/90°
- Leg up and down 20°/100°
- Back up and down 90°/40°
- Trendelenburg and reverse 30°/30°
- Power kidney bridge height 6 inches, (6700B is 5-3/4 inches high).
- Secondary controls on base of table
- Return to Level (Anatomically correct return to horizontal position from articulated positions)
- Battery life indicator
- Removable leg and back sections
- Table weight: 900lbs.



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.