



auxoMEDICAL

Skytron Infinity Series- Fully Refurbished

The Skytron® Infinity® Surgical Light is a powerful light suitable for general procedures, both in a hospital or small operating room setting. The light provides high quality shadow control and superior focusing capabilities with both a sterile central handle control and non-sterile side controls. The glare reducing soft white light is ideal for clarity and a pristine viewing field, crucial to a successful procedure.

Features:

- Electrical: 120v, 60 Hz, 1 phase
- Power consumption: 570 to 600 watts (1,950 to 2,448 BTU/hr)
- Amperage: Approximately 5.0 amps, fuse protected
- Color temperature: 4,000°K
- Beam temperature: Less than 25,000 microwatts per square centimeter as recommended by the I.E.S. handbook
- Bulb type: Halogen 24v/50w (operates at 18.5 to 20v)
- Weight: 260 lbs
- Moment load: 640 ft lbs
- Ceiling height: Min 9'4", Max 10'3.5"
- 22" in diameter
- Multi-lamp optics with True Focus Capability™ for greater working distances
- Glare reducing soft white light
- Self-contained wall-mounted intensity controller
- Depth of Field: 20" .
- Sterilizable handles are optional
- Center Mounted Suspension System
- True Focus Capability For Greater Working Distances
- Variable Intensity Controller



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.