



auxoMEDICAL

Statim 2000 Sterilizer - Fully Refurbished

The STATIM 2000 is a reliable cassette sterilizer that's compact enough to fit in procedure areas and fast enough to sterilize jewelry, expensive Lasik, PRK, endoscopic and cataract instruments, including phaco hand pieces, between procedures – it can complete an unwrapped sterilization cycle in just six minutes and a wrapped cycle in 14 minutes. The STATIM 2000 Classic boasts a 6-minute cycle, but has the added advantage of employing a fully removable chamber so you can aseptically transport a set of freshly sterilized instruments directly to the procedure area.

Features:

- No plumbing required. Integrated reservoir.
- Front panel controls.
- Low water cut-off. This prevents burnout of boiler elements.
- Safety switch automatically activates when pressurized
- Chamber Volume: 1/2 Gallon (1.9 Liters)
- Reservoir: 1 Gallon (3.8 Liters)
- Cycle Time: 6 minutes
- Fully automatic
- Three sterilization programs
- Advanced microprocessor controls all functions and maintains optimal sterilization



Specs:

- External Size (LgthxWdthxHgth): 17 3/4" x 16 1/4" x 6"
- Cassette Internal Dimensions: 11" x 7" x 1 1/2"
- Reservoir Capacity: 4 Liters
- Weight Without Water: 46 Pounds
- Power Consumption: 1300 Watts
- Cycle Time:
 - Unwrapped: 6 Minutes*
 - Wrapped: 12-15 Minutes*
 - Rubber & Plastics: 20 Minutes*

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.