



auxoMEDICAL

Statim 5000 Sterilizer - Fully Refurbished

The SciCan Statim 5000 cassette autoclave (or steam sterilizer) uses a steam injector system provides temperature consistency and expels the air which causes oxidation. The Statim 5000 has an additional sterilization cycle and a roomier cassette designed to process bulkier instruments and greater loads. It's a reliable cassette sterilizer that's compact enough to fit in patient treatment areas and fast enough to sterilize expensive Lasik, PRK, endoscopic and cataract instruments, including phaco handpieces, between procedures. SciCan's Statim cassette sterilizers have a biological effectiveness that has been proven in tests for hollow and solid instruments at many internationally recognized institutes in Canada, the United States and Europe.

Features:

- Wrapped Sterilization cycle: 17.5 min
- Unwrapped Sterilization cycle: 9 min
- Fully Automatic
- Four sterilization programs
- Advanced microprocessor
- Fast, Compact, Reliable
- Effective Every Time
- Gentle on Instruments
- Easy to Use and Maintain
- Exceptional instrument turnover
- 15 minute Dri-TecR Drying Technology
- Prolongs instrument life
- Reduced instrument investment



Specs:

- Unit Dimensions: L 55.5 cm x W 41.5cm x H 19cm | L 21.75" x W 16.25" x H 7.5"
- Weight Without Water: 30kg / 66lbs
- Cassette Internal Dimensions L 38cm x W 18 cm x H 7.5 cm | L 15" x W 7" x H 3"
- Power Consumption: 110-120V 60Hz, 11A, 1300 W
- Reservoir capacity: 4.0 litres (distilled water)

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.