



## Steelco DS 500 - Fully Refurbished

The Steelco DS 500 washer/disinfector is ideal for doctors' offices, surgery centers and small hospitals to wash, disinfect, and dry surgical instruments. Auxo Medical can easily install this unit under (or on) a counter as single door version. The high capacity chamber the washing and hot air drying on two independent levels, which are normally seen on larger machines, allow this device to be used following on everyday needs. The powerful cleaning and disinfecting system guarantees short programs duration that reach high level performances in terms of cleaning quality, thermodisinfection and drying.

### Features:

- Heavy washing pump to ensure high water flow rate combined with effective spray pressure.
- Filtered forced air drying system with adjustable time and temperature setting, helps to ensure a complete drying, inside and outside, of all the surgical instruments and tubes.
- Washing and DI disinfection temperature are fully adjustable up to 93°C. Temperature is monitored by two independent sensors.
- Standard steam condenser prevents vapours from entering into the washing area.
- Inner cabinet, washing arms and tank filters made of high quality AISI 316 L stainless steel.
- Two automatic liquid dispensers.
- Triple water filtering system captures residue, preventing re-circulation and extending the pump life



### Specs:

- Dimensions WxDxH: 23.62" x 24.80" x 33.46"
- Chamber volume ~171 lt / 45.17 Gal US
- Basket volume ~151 lt / 39.89 Gal US

### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.