



Steris 4085 - Fully Refurbished

The STERIS 4085 General Surgical Table's design gives a surgical team the flexibility, control and convenience they need for today's patients and procedures. The Steris 4085 general surgical table handles virtually all the surgical specialties you expect a general surgical table to handle and provides outstanding access for C-arm, permitting high quality images crucial to cardiothoracic, orthopedic and neurosurgical procedures.

Features:

- Four section tabletop
- Self-leveling floor locks
- Removable leg section
- Powered, radiolucent kidney elevator
- Easy to use, easy to read hand control with digital data window
- Prevents sectional and table base collisions with the Auto Limit Sensor™

Specs:

- Dimensions: 81"(L) x 20"(W) x 26" - 45"(H)
- Table Slide Range 9" (227 mm) to head, 9" to foot (18" total)
- Patient Weight Capacity 1,100 lbs. (500 kg)
- Table Shipping Weight 560 lbs. (254 kg)
- Trendelenburg / Reverse 30° / 30°
- Lateral Tilt (left/right) 20° / 20°
- Head Section +90° to -90°
- Back Section +80° to -40°
- Leg Section 0°/-105° (removable)
- Flex/Reflex 140° / 100°
- Manual Override Yes
- Perneal Cut-out Yes
- Kidney Elevator 4" powered radiolucent



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.