



auxoMEDICAL

Steris ASC 2000 OR Table - Fully Refurbished

The ASC 2000 Ambulatory Surgical Care Table is STERIS's premier surgical table for the Ambulatory Market. This table is designed to provide complete and flexible positioning for outpatient procedures. The unique, offset table column design facilitates surgical team access and comfort while providing superior C-Arm access.

Features:

- Latex-free tabletop pad
- Bakelite tabletop
- Power adjustable height, back (upper/lower), Trendelenburg and lateral tilt functions
- Manual, removable, split leg and headrest sections
- Pendant hand control
- Mobile base with manually operated floor locks
- 400 pound (181 kg) patient weight capacity
- Up to two weeks operation on internal battery power

Specifications:

- Dimensions: 22" x 79" x 27.5" to 41"
- Height: 27.5" to 41"
- Weight: 485 lbs.
- Trendelenburg Range: 45° ±1° Reverse Trendelenburg Range: 20° ±1°
- Lateral Tilt Range: 20° ±1°
- Lower Back Motion: UP: 75° ±1° / DOWN: 15° ±1°
- Upper Back Motion: UP: 40° ±1° / DOWN: 55° ±1°
- Head Motion (Manual): UP: 33° ±1° / DOWN: 87° ±1°
- Leg Motion (Manual): UP: 0° / DOWN: 90° ±1°

Additional Options: Armbands and pads, Straps, Positioners, Urology Extension, Narrow Headrest, Stirrups, Prepping Positioners, Urology Drainage bags, Hand Tables and Traction Units



Why an Auxo Medical refurbished Surgical Table?

Auxo Medical tables are refurbished to the highest standards in the industry. We disassemble each table, make necessary cosmetic upgrades, replace hydraulic components and seals and return the table to original OEM Specifications. We do not substitute cheaper non-OEM parts like many of our competitors. Each table is reassembled, tested and is warranted against defects.