



auxoMEDICAL

Steris Amsco 3080- Fully Refurbished

The AMSCO 3085 General Surgical Table is one of the most reliable and popular surgical tables. Trusted for years by healthcare facilities across the world, the 3085 is a cost-effective solution for the general OR. The Amsco® 3080 Surgical Table features a radiolucent top and can be adjusted to the same variety of positions as the 2080 model. It comes standard with a remote control for easy and efficient patient positioning. This is one of the most durable, flexible, and reliable surgical tables available.

Features:

- Remote control
- C-ARM Compatible
- Return to level and special positioning models available
- Positions include: Trendelenburg, back (raise or lower), side (lateral tilt), foot (raise or lower), flex ("v" or inverted "v")
- Up to 500 lb capacity
- Durable, flexible and reliable
- Brings more efficiency to the operating room
- Fluoroscopic capabilities
- Simple patient positioning
- handles urology and neurology, pediatrics, ENT and endovascular, and even orthopedics, including skeletal traction, with accessories.
- Mobile, electro hydraulically operated
- Floor Locks
- Length: 76 in (1930 mm) Width: 20 in (508 mm)
- Weight of operating table: 920 lbs (417 kg)
- Height Adjustment: 27 - 44 in (686 - 1118 mm)
- Lateral tilt (left/right): 18°
- Trendelenburg/Reverse Trend: 25°/25°
- Back up/down: 55°/25°
- Leg up/down: 80°/105°



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.