



## Steris Amsco Century Small Series- Fully Refurbished

The Century® Sterilizer meets the needs of today's O.R. and other areas of the hospital. The Century® features foot pedal operation of the automatic sliding door for hands-free loading and unloading. The vacuum fluorescent tough panel permits direct operator cycle selection. The large display of time remaining during the cycle allows the operator to perform other tasks while awaiting cycle completion.

### Features:

- Chamber material and pipe work 316L stainless steel
- Fully jacketed chamber
- Vertical opening doors with glass panel
- ALUS chamber rail for racks
- 6,5" (165 mm) high resolution color touch screen
- Interface to tracking system XLM and PDF format
- External USB connection for data download (on loading side)
- Matrix printer on loading side with A4 print out capability
- Additional touch screen provided on the unloading side
- Independent temperature probe (chamber)
- Pneumatic valves
- Air filter
- Pressure gauges for jacket, chamber, steam generator (loading side), chamber (unloading side)
- Water level gauge for steam generator
- Safety temperature switch for door opening
- Acoustic signal for end cycle and alarms
- Programmable timer start for test cycles
- Potential free outputs for critical alarms
- Sanitary piping with Tri-clamp
- Water Saving System For Vacuum Pump
- Separate water feed line for steam generator and vacuum pump
- Drain discharge cool down



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.