



Steris Amsco QDJ Series - Fully Refurbished

Steris Amsco QDJ Series warming cabinets are designed to store and warm sterile intravenous (IV) and surgical irrigation solutions, linens, and blankets to programmed temperature settings for patient applications. The digital controls provide an LED readout that verifies set-point and actual unit temperature for each compartment. Temperatures are maintained at minimum cost with the addition of a visual indicator for when the door is ajar. A heating element interruption when the door is ajar prevents accidental overheating of supplies and saves on utility costs.

Features:

- LED readout for instant temperature confirmation
- Superior durable materials to withstand years of use
- Adjustable shelving for flexible organization and storage
- A design to meet AORN Recommended Practices for Safe Environment of Care, including technology to warm IV fluids (RP VIII.b.3)
- Reliable heating with uniform and efficient heat transfer
- Visual access to inventory with dual-pane tempered glass or traditional solid steel doors
- Temperature lock-out for consistency and control
- Celsius and Fahrenheit selection for convenience and accuracy
- Door Option: Steel or dual-pane tempered glass
- Installation: Freestanding or recessed
- Temperature Range: 90°F (32°C) to 160°F (71°C)



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.