



auxoMEDICAL

Steris Amsco V136 - Fully Refurbished

Completely refurbished, 24" x 36" x 36", single door sterilizer recess mounted, or optional cabinet package (G136 variant for gravity operation) with operating temperature of 250-270F. Designed for fast, efficient sterilization of heat and water-resistant materials with the same capabilities as a gravity sterilizer. The V136 pre-vacuum sterilizer is equipped with pre-vac, gravity, flash, express, leak test, and daily air removal test cycles. The Amsco V136 sterilizer comes complete with a one year parts and labor warranty (6 months on preventive maintenance part), five year warranty on the pressure vessel.

Features:

- Florescent display for easy reading
- Ink on paper impact printer
- 24" x 36" x 36" chamber
- Rack and 2 shelves
- Recessed
- Non-lubricated, steam-activated door seal
- Help screens for programming/troubleshooting alarm conditions
- Electronic water conservation valve
- Pre-vacuum Steam Sterilizer
- Emergency manual exhaust valve
- Communication interface compatibility w/ most PC devices



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.