



Steris Manual Dual Bay Flexmatic Scrub Sink - Fully Refurbished

Dual Bay Flexmatic scrub station without timer or water control sensor. Designed in accordance with Underwriters Laboratories safety standards, AMSCO Flexmatic Scrub Sinks fulfill the requirements of any operating room. Available with 1, 2, or 3 stations, these sinks offer hands-free operation of both water and soap. Safety features minimize contamination risk, and design options and accessories provide convenience to your staff. . Each unit is constructed of welded 16 gauge stainless steel with a matte finish for easy cleaning and low maintenance. The built-in thermostatically controlled mixing valve provides optimum water temperatures up to 115° F (46° C). These durable sinks feature continuous welded compartments and coved corners constructed entirely of stainless steel for easy cleaning. Easy-to-use knee panel makes water control quick and convenient.

Features:

- Two Knee Push Pedals for Activation
- Chair Carrier Included
- Installation not included
- Soap Dispensers not included
- Dimensions: 64 x 27 9/3 x 5 1/2"
- Stainless steel construction with welded compartments and coved corners for easy cleaning
- Standard "rose" spray head produces a shower head type pattern without aerating the water
- Hands-free soap – 2 cc. of soap are released with each push of the foot pump or knee operated dispenser



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.