



Steris SurgiGraphic 6000 - Fully Refurbished

The SurgiGraphic® 6000 Image Guided Surgical Table offers a new level of operating room performance and control – allowing radiology-quality imaging in the operating room. The table delivers the patient positioning and mobility you demand for most minimally invasive procedures. The SurgiGraphic surgical table enables radiology-quality imaging for minimally-invasive vascular, orthopedic, or other procedures requiring fluoroscopy; lateral and longitudinal float, Trendelenburg and reverse Trendelenburg, lateral tilt, height adjustment, and the most useful combinations of these movements.

Features:

- Physician controlled Manual 4-Way Float Palm Handle for immediate and exact patient positioning without moving the C-arm.
- Simultaneous Lateral Float and Trendelenburg/ Reverse Trendelenburg.
- Simultaneous Longitudinal Float and Lateral Tilt.
- Height adjustable
- Lateral Tilt: 16°
- Trendelenburg/ Reverse Trendelenburg Tilt: 16°
- Customized hand control with programmable memory positions.
- 500 lb. patient weight capacity
- Patent pending specially designed detachable accessory side rail clamp system for complete metal-free imaging after removal. Designed to hold retractors, arm boards and all O.R. accessories.
- Concealed swivel casters allow for superior maneuverability.
- Self compensating floor locks
- Power plugs into either front or rear of table
- Stainless steel base cover
- UL listed
- Dimensions: 92"(L) x 21"(W) x 32" - 42"(H)



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.