

## Stryker 1061 OB/GYN Stretcher - Fully Refurbished

---

The Stryker 1061 can go from a standard transport stretcher to a stretcher that fits every OBGYN need. Go straight from performing pelvic and OB exams to general use with the quick fold hideaway stirrups and collapsible foot extension. Don't buy multiple stretchers for different tasks. Own one stretcher that does it all.

### Features:

- 4" (10 cm) Ultra Comfort Mattress with sealed seams
- 8" (20 cm) Omni Surface casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics with constant rate descent
- Dual side-mounted foot control hydraulics + retractable fifth-wheel steering
- Four-wheel, steel-ring brakes with dual side and head activators
- Integrated storage tray with oxygen bottle holder
- Permanently-attached folding foot section
- Pop-up push handles, roller bumpers, Glideaway side rails

### Specifications:

- Overall Length 81.5" (207 cm)
- Overall Width
  - Side rails Up 31.5" (80 cm)
  - Side rails Down 26.75" (68 cm)
- Weight Capacity 500 lbs. (228 kg)
- Height Range (to litter top)
  - High 36" (91 cm)
  - Low 22.5" (57 cm)
- Patient Surface 26" x 77" (66 x 196 cm)
- Side rails 13" x 39" (34x99cm)
- Caster Diameter 8"



Accessories: Calf Supports, Brackets, Defibrillator Tray, IV Caddy and Poles, Sidereal Pads, Patient Restraints, Serving/Instrument Tray

### Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rebuild or replace all worn parts or cylinders, and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.