



Stryker 1080 Fluro Stretcher - Fully Refurbished

Access, functionality, and safety are the hallmarks of Stryker's Fluoroscopy Stretcher. The specialty platform integrates multiple advanced features for enhanced application and control. The stretcher's features are designed to simplify operation and minimize the need to re-position the patient during fluoroscopic imaging. The full-length C-arm access, fold down siderails, and dual foot controls increase flexibility and ease of imaging.

Features:

- Weight Capacity: 500 lbs
- Radiolucent Litter Top
- Dual Fibresin Pneumatic Backrest
- Dual End Hydraulic Foot Controls
- 6" Non-Conductive Casters
- Fold Down Side-Rails
- Four Wheel, Steel Ring Brake System with Dual End Activators
- Four IV Receptacles
- 3" Thick Stretcher Pad
- Completely Patient Ready

Specs:

- Overall Bed Length: 87"
- Overall Bed Width: 30.5"
- Manual Height Adjustment: 24" - 34"
- Patient Surface Dimensions: 25.5" x 84"
- Side-Rail Dimensions: 13" x 55"
- Fluoroscopy Window: 19.5" x 66" (50 cm x 168 cm).
- C-Arm Access: 54" (Length) (137.16 cm).



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.