



Stryker 5050 Stretcher Chair - Fully Refurbished

The Stryker 5050 Stretcher Chair offers multiple position options for superior versatility. The chair's easy positioning from horizontal to chair position helps achieve superior clinical performance and patient comfort. Ideal for surgery centers and small spaces, the 5050 maximizes patient throughput.

Features:

- Adjustable height footrest
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Fixed position pushbar (model 5050 only)
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual side activators
- Hands-free Trend.
- Swing-down siderails with comfort armrest
- Thermoformed ABS hood with integrated storage tray and oxygen bottle holder

Specs:

- 2" New vinyl mattress
- 6" casters
- 24" wide patient surface
- Height Range: 33.5" high 22" Low
- Patient Surface: 24"x74"
- Overall Width: 30"
- Litter Positioning: Backrest 0-90°
- Leg 0-80°
- Trend/reverse Trend: +/- 18°



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.