



## Stryker 735 Transport Stretcher - Fully Refurbished

The simple design of the Stryker 735 Transport Stretcher provides safe and efficient transport for an economical value. Ideally suited for the ASC environment, this stretcher features fifth wheel mobility for easier traction and cornering. Additionally, brakes on all four casters lock both swivel and rotation for superior security.

### Features:

- 3" (8 cm) mattress
- 8" (20 cm) Omni Surface casters
- 26" (66 cm) wide patient surface
- Fold-down siderails (fold to head)
- Four IV receptacles
- Four-wheel, steel-ring brake system with dual end activators
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Roller bumpers
- Stationary foot-end litter
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side-mounted foot control hydraulics
- Hands-free Trendelenburg
- Overall length: 83"
- Overall width: 30"
- Weight Capacity: 500 lbs.



### Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.