



Tuttnauer 2340EA EZ9 Autoclave - Fully Refurbished

This table-top autoclave is designed for the sterilization of wrapped and unwrapped instruments, glassware and other related items found in dental, medical and veterinary clinics, tattoo parlors, first aid rooms, hospitals, laboratories and other related industries. The enhanced closed door active drying system maintains the sterile field and ensures complete drying of packs and pouches. The efficient air pump forces hot air circulation after the sterilization cycle for faster, more efficient closed door drying

Features:

- Front panel controls
- Safety device to prevent premature opening while unit is pressurized
- Time saving Reset button
- Dual safety thermostat to protect against over heating
- Automatic shut off when sterilization and dry cycle is completed
- Electronic Digital Control
- Front side drain valve for ease of releasing reservoir water.
- 5 Gallon Chamber (19 Liters)
- Holds 3 trays
- Cycle Time: Varies based on items to be sterilized - Wrapped instruments 20-30 minutes



Specs:

Temperature Range: 250 degrees F thru 274 degrees F (121C - 134C)

Electrical: 120 volts, 1400 watts, 50/60 HZ, 11.7 AMPS

Unit Weight: 90 lbs (40.8kg)

Overall Size:

20"(W) x 21-1/2"(L) 14-1/2" (H)

Chamber Diameter:9" (23 cm)

Chamber Length: 18-1/2" (46.9 cm)

Tray Dimensions: 6-3/4"(W) x 16-1/4"(L) x 3/4"(H)

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.