



Tuttnauer 3870 - Fully Refurbished

This is Tuttnauer's largest manual table-top autoclave. Believe it or not there is still a demand for the 3870m analog unit among many healthcare facilities and veterinary practices. This beefy, no frills autoclave covers all the basics needed for proper sterilization. It can always be counted on even under harsh conditions. The large 15" x 30" chamber accommodates longer instruments or a large number of smaller instruments, making it a perfect match for spay and neuter clinics, equine practices, busy animal hospitals and multi-practice medical offices.

Features:

- A long life electro polished chamber and door
- Double safety locking device prevents door from opening while chamber is pressurized
- Overheat cutoff prevents burnout of heating elements in case of insufficient water in the chamber
- Drain valve is located on the front
- Dual safety thermostat to protect against overheating
- International certifications
- Closed Door Dry Cycle
- Freestanding
- Filled with distilled water



Specs:

- Chamber Dimensions: 15" x 30"
- Volume: 22 gal
- Overall Dimensions D x W x H: 34.5" x 26" x 20.6"
- Standard Cassettes: 15 full
- Standard Unwrapped Cycle Time
 - Cold: 55 mins
 - Hot: 40 mins
- Tray Dimensions D x W x H:
 - 11.3" x 26.6" x 1"
 - 13.8" x 26.6" x 1"

Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.