



UMF Bariatric Exam Table - Fully Refurbished

The 5100 series line of precision built exam tables come standard with the most innovative and requested feature in the industry. UMF's high quality doesn't come at the expense of practicality or affordability. Every UMF exam table is designed to provide maximum patient comfort and years of trouble free use. There's a reason the 5100 series is the most recognized and imitated line of exam tables in the industry.

Features:

- UMF Ultra Comfort Top (vacuum formed and seamless)
- Easy-glide comfort stirrups (3 positions)
- Upholstered removable leg rest
- Seamless table top
- Pull-out all steel non-slip step (7" high)
- Paper roll holder
- Plastic 9" x 9" drip pan
- Hospital grade receptacle
- Two front storage drawers
- Three side storage drawers
- 500 lbs capacity

Specs:

- Overall: 27" W x 58.5" – 77.5" L x 31.75" H
- Stirrup: 3 lateral positions 22 degree, 32 degree, 42 degree, infinitely adjustable over length of the bar
- Front drawers: (2) 19" W x 16" D x 4.25" H
- Side drawers: (3) 19" W x 16" D x 4.25" H
- Pelvic lift 7.5 degrees



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.