



## Walkmed 350VL Infusion System - Fully Refurbished

The WalkMed 350VL pump is indicated for intravenous, subcutaneous, arterial, enteral, and epidural infusion of: antibiotics, analgesics, chemotherapeutic agents, and other medications or fluids requiring precisely-controlled infusion rates. The delivery system in WalkMed ambulatory infusion pumps is highly regarded for its trouble-free performance and its ability to deliver accurately, even at low infusion rates. The WalkMed 350VL is ideal for chemotherapy.

### Contraindications for Use:

- Infusion of blood and blood products
- Infusion of insulin
- Infusion of critical medications whose stoppage or interruption would cause serious injury or death
- Use in ambulatory regimens by patients who do not possess the mental, physical, or emotional capability to operate the pump properly: or who are not under the care of a responsible individual

### Features:

- Patient Lockout
- Volume Tracking
- AutoOcclusion Alarm Recovery
- Symbolic Keypad
- Dimensions: 4.4 x 4.0 x 1.5 in
- Weight: 12.8 oz. / 360g

### Safety & Convenience:

- Keypad Lock-out for Patient Protection
- Auto-Occlusion Alarm Release for Sensing & Silencing Alarms
- Programmable Total Volume Limit & Delivery Tracking



## Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.