



auxoMEDICAL

Zimmer ATS 3000 - Fully Refurbished

The A.T.S.® 3000 Automatic Tourniquet System is a dual-port, dual-cuff system with microprocessor controls and dedicated ports for supplying and measuring pressure independently. With the innovative Limb Occlusion Pressure (LOP) feature, the A.T.S.® 3000 Tourniquet combines the latest in advanced tourniquet technology with the well-established Zimmer® tradition of safety, reliability, and convenience.

Features:

- Dual Port, Dual Cuff System
- Limb Occlusion Pressure feature that detects, calculates, and reports cuff pressure necessary to achieve complete blood occlusion in the operative limb.
- Four Hour battery back up that maintains pressure if power is lost
- Visual and Audible Alarms
- Self check calibration
- Lightweight and portable
- Ports for supplying and measuring pressure independently
- Cuff Lockout feature that tells user to confirm the deflation command during procedures to reduce complications due to sudden cuff deflation
- Color-Coded Cuffs: From 8 - 42"
- Carrying Handle
- Independent pressure monitoring
- Rapid inflation and deflation
- Built in safety features
- Positive locking connectors for leak free inflation
- Lightweight unit is easily used on a table or IV pole



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.