



# auxoMEDICAL

## Zoll E Series - Fully Refurbished

ZOLL Medical's full-featured defibrillator was designed specifically for the rigors of the EMS environment. No other defibrillator offers the configurable options of the E Series. No other defibrillator is as rugged and ready. The E Series offers ZOLL's Real CPR Help®, which measures chest compressions and rate and depth in real time, and provides visual and optional audible feedback. All CPR data can be recorded and reviewed using RescueNet Code Review software. Documentation is made easy with data transfer options that take the "paper" away. The See-Thru CPR® allows CPR artifact to be filtered, letting you see organized rhythms without pausing compressions.

### Features:

- Patient Connection: 3-lead, 5-lead or 12-lead ECG cable, paddles or MFE pads, selectable by front panel switch.
- Measures chest compressions and rate and depth in real time, and provides visual and optional audible feedback
- Unique feature allows CPR artifact to be filtered, letting you see organized rhythms without pausing compressions.
- New innovative roll cage design and shape.
- Screen usable in pitch darkness or direct sunlight.
- Easily and accurately diagnoses CO poisoning with the push of a button
- Built-in GPS clock, allowing users to synchronize dispatch, defibrillator, and intervention call times, improving overall data accuracy.
- Compatibility with ZOLL's sealed lead acid 4410 Smart Battery as well as rechargeable lithium ion.
- ZOLL's proprietary Rectilinear Biphasic waveform
- Dimensions: 5.75 in. high x 13.1 in. wide x 10.5 in. deep
- Weight: Approximately 13.2 lbs (5.99 kg)



### Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.