



auxoMEDICAL

Zoll M Series - Fully Refurbished

The Zoll® M Series™ is half the weight and size of any unit in its class. Provides everything you need for cardiac resuscitation including easy-to-use defibrillation, external pacing, and advanced monitoring capabilities. Plus, it gives you these important options: pulse oximetry and capnography, a sophisticated code marker and data collection system, and summary report with printer.

Features:

- The smallest, lightest, brightest defibrillator unit in its class
- Easy handling and optimal viewing
- Built-in AC power and battery charger
- Consistent, straightforward controls
- Intuitive menus lessen training costs, minimize operator confusion, and reduce the chance for error
- ZOLL's proprietary Rectilinear Biphasic™ waveform has lower energy and current
- 40-msec Rectilinear Pacing waveform provides superior efficacy while optimizing patient safety
- Patented external Pacing offers superior capture at lower mean thresholds
- Offers defibrillator event documentation via internal memory
- Removable PCMCIA cards for easy data transfer and archiving
- M Series communication capability makes record keeping a breeze
- Pre-connection capability and common connectors help reduce stress in critical situations
- 12-Lead Interpretive ECG
- Non-invasive Pacing
- Multifunction Cables or Integrated Hard Paddles
- Size: 6.8 in (17.3 cm) high x 10.3 in (26.2 cm) wide x 8.2 in (20.8 cm) deep.
- Weight: 11.5 lb. (5.23 kg) with Multi-Function Cable and battery; 13.5 lb. (6.14 kg) with paddles.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.