



BK Medical Falcon 2101 EXL - Fully Refurbished

The Falcon 2101 EXL's high image quality and versatility makes it useful for a variety of applications. With features like intuitive interface and user-defined settings, the Falcon optimizes your imaging conditions as well as the use of your time. The sleek, compact and mobile design provides ultrasound where you need it, precisely when you need it. And, it is equipped with the right tools: a comprehensive collection of transducers with corresponding puncture and biopsy equipment.

Features:

- Frequency range of 1-20 MHz provides the highest resolution in its class.
- Four transducer connectors (three array and one mechanical).
- B-K Medical's patented Continuous Uniform Focusing technology ensures optimal resolution at all image depths.
- Compact and mobile, so it can be used wherever and whenever you need it.
- Ergonomic, with height-adjustable keyboard and monitor and logically grouped controls.
- Full 15" flicker-free monitor provides optimum viewing.
- Optional fully integrated 3-D package, including tracked pullback support for relevant transducers.
- Simultaneous split screen feature lets you display two active ultrasound images on the monitor at the same time.
- Convenient storage space on the unit itself, with a compartment under the keyboard for accessories.
- Upgradability ensures you of long-tem investment security.



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.