



Brewer Element Treatment Table - Fully Refurbished

The Brewer Element Treatment Table offers fundamental features every practice needs. A robust steel frame construction ensures maximum stability and safety for patients while its contemporary design adds style to any healthcare facility. Ideal for basic examinations, pediatrics, orthopedics, physical therapy, industrial areas, and school nurse stations.

Features:

- Solid, all steel frame for durability
- Two piece design for freight savings
- Contoured, seamless upholstery for easy cleaning
- Variety of configuration options available
- 500 lbs weight capacity
- All steel frame construction with chip resistant powder coated finishing
- Standard concealed paper roll holder accepts 21" & 18" rolls (53.3 cm & 45.7 cm rolls)
- Optional drawer models can be made left or right handed with space for up to six rolls
- Adjustable back models have six positions
- Two-piece top design ensures easy assembly, single person delivery, and freight savings
- Contoured, seamless upholstery for easy cleaning
- 2" (5.1 cm) thick foam for patient comfort
- Height: 31" (78.7 cm)
- Length: 74" (187.9 cm)
- Width: 28" (71.1 cm)



Why an Auxo Medical Refurbished Table or Stretcher?

Auxo tables and stretchers are refurbished to the some of the highest standards in the industry.

We completely disassemble each table and stretcher, and completely rework or replace all worn parts and return each table and stretcher to OEM standards. Each table and stretcher is reassembled, tested and is warranted against defects.