



## Multi-Use Imaging Table with Stirrups and Drop Window

---

Clinton offers a broad range of medical tables designed to meet the needs and preferences of physicians, hospitals and facility managers. From classic sturdy wood treatment tables to steel and cabinet models, Clinton tables are built for ease of use, durable performance and patient comfort.

### Features:

- Three-section table with pneumatic adjustable backrest, manual footrest, stirrups and window drop
- Footrest is adjustable from either side and raises up to 30° from flat position
- Optional side rails & casters (see opt
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Power supplies and motors meet UL 60601-1 and CAN/CSA C22.2 601.1
- Industry leading 5-year, limited warranty
- 34" wide model available
- 550 lbs load capacity under normal use

### Options:

- Power Backrest
- Headrest/Power Tables
- Safety Straps
- 2-Hook IV Pole & Mounting Bracket
- Locking Casters
- Power Footrest
- Safety Rails



### Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.