



Drager Apollo Anesthesia Machine- Fully Refurbished

The role Apollo plays in optimizing your anesthesia workflow is difficult to understate. This medical equipment decreases process costs with features like the fully automatic self test. Since user interaction is not required, practitioners become freed up to take on more tasks. In addition to providing best-in-class therapy, Apollo also gives you an open-platform architecture – enabling this machine to be readily integrated into your overall workflow.

Features:

- User friendly ergonomics
- Automated startup
- Advanced ventilation capability provides ICU quality ventilation in the OR
- Spontaneous breathing support
- Electrically powered piston ventilator requires no drive gas
- Enhanced monitoring capability provides user with the information they need to make clinical decisions

Specifications:

- VT 20-1400 mL
- Frequency: 3-80 bpm
- I:E ratio 5:1 (max)
- Inspiratory Pause 0-60%
- PEEP to 20 cm H₂O
- Adjustable Peak flow, plateau time
- Pressure limit: 70 cm H₂O
- Inspiratory flow (pressure control mode) up to 150 L/min

Additional Options: D-Vapor / Vapor 2000, Infinity Omega Solution, Innovian Anesthesia



Why purchase Auxo Medical Reconditioned Biomedical Equipment?

Auxo thoroughly inspects and re-certifies all of its reconditioned biomedical equipment. Batteries are checked or replaced, and each unit is cosmetically repaired. Calibrations are made and all safety inspections will be completed. In most cases our equipment warranties meets or exceeds those of the OEM.