



Hausted All Purpose Chair Hydraulic Models - Fully Refurbished

The Hausted All Purpose Chair (APC) is a multi-position chair and stretcher intended for intra- and inter- departmental transport of patients within a health care facility. The APC can be positioned to allow for patient treatment, transport or recovery.

Features:

- Caster diameter - 5" (127 mm)
- Fold-down foot rest (height adjustable)
- Padded armrest:
 - Height - 12" (305 mm)
 - Length - 29" (737mm)
- Polymer-coated support surface
- Four universal accessory wells
- Floor level foot pedals for height adjustment
- Backrest and Trendelenburg release levers
- Degree of movement: Backrest - 0-80° and Trendelenburg - 10°
- 3" (76 mm) comfort cushions
- Built-in chart holder
- Built-in push bar
- O2 Tank holder
- Restraint/safety straps
- Oxygen tank holder
- Removable IV rod
- Weight Capacity: 325 lbs.



Why an Auxo Medical Refurbished Bed or Stretcher?

Auxo beds and stretchers are refurbished to the some of the highest standards in the industry. We completely disassemble each bed and stretcher, and completely rework or replace all worn parts and return each bed and stretcher to OEM standards. Each bed and stretcher is reassembled, tested and is warranted against defects.